Release Notes

coDiagnostiX 10.1

English

New surgical guides

While earlier versions of coDiagnostiX were solely intended for dental implant planning, beginning with version 10, coDiagnostiX is intended for planning in dentistry. This opens the software for a series of new applications, while the main workflow continues to be the same.

Bone reduction guides

Since in certain cases, bone reduction is necessary prior to dental implant placement, coDiagnostiX 10 now offers the planning and design of bone reduction guides as part of a collection of new features for the edentulous workflow. These guides support cutting of osseous structures by guiding an instrument (e.g. piezo saw) along a predefined profile.

coDiagnostiX supports the definition of a cut profile. This profile then serves as a basis in the regular guide design wizard to design a bone reduction guide.

This procedure can also be applied to external sinus lift, apicoectomy, and gingivectomy planning.

Endodontic guides

With endodontic guides, we now offer a smart and safe treatment approach for cases with pulp canal calcification or root canals that are difficult to locate. Endodontic guides accurately guide the drill to the obliterated root canal to create an access cavity.

For this, the coDiagnostiX database now includes endodontic drills. Endodontic planning in coDiagnostiX is realized by selecting endodontic drills from the integrated database and placing them in the same way as implants. Subsequently, appropriate sleeves are assigned and a surgical guide designed with the regular coDiagnostiX guide design wizard.

Model scan import

We improved the handling of imported surface scans and other 3D objects which are not optimal for treatment planning to better aid virtual tooth extraction and mesh optimization.

This new functionality can be found in the Model Scan Import dialog and in the context menu of imported model scans and other 3D objects.

Other

- The loading speed of the implant selection dialog has been significantly improved.
- When a lateral pin is added, it is now automatically placed at the selected tooth and at an angle. This ensures that the pin is in the best position to start fine alignment. In addition, pins are now assigned an undefined tooth position (XX) to prevent them from blocking the tooth position for implants.
- A new printout has been added to allow users more customization options. It contains all images stored in the Image Management for a particular patient.

Improvements in EASY mode

- It's now possible to measure the distances in the Implant Alignment step in EASY mode. Previously, this function was available only in EXPERT mode.
- The Dataset Preparation step in EASY mode now features a segmentation threshold slider. It allows users to optimize the appearance of their 3D reconstruction from scan data more easily.

Bugfixing

- Increased stability of the surgical guide designer.
- Fixed a rare problem where producers were not able to read datasets from their locked clients.
- Fixed a rare problem that occurred in the Straumann surgical protocol print-out. Sometimes, the order of tooth positions did not correctly match the dental notation system that was chosen.
- Fixed a problem with five CAMLOG implants, where the visible diameter of the 3D object was 0.4 mm smaller than the value from the catalog. The following implants were affected: K1054.4311, K1054.4313, K1054.4316, K1053.4316, and C1086.3809.

Known issues

- coPeriodontiX is not accessible until further notice.

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