Release Notes

coDiagnostiX 9.12.1

English

Virtual planning export improved

When converting scans into, for instance, models with implant analogs, the virtual planning export now allows users to store the resulting object directly in the current case planning. Users no longer have to re-import an STL file. This new and convenient import option applies to intraoral, model or impression scans as well as to segmentations and all other user-imported 3D objects.

Support for digital marker trays

coDiagnostiX now supports the digital marker trays for the DENACAM system from Mininavident.

Refined implant-abutment alignment

An issue regarding abutment alignment has been fixed for implants with a 30% angled shoulder.

Rapid Shape RSHAPESTL output file format added

Printing guides on a Rapid Shape printer has become easier. When exporting a surgical guide to the printer, the RSHAPESTL file format now transfers meta data, which helps accelerate the printing process.

Surgical protocol improvement

We improved the surgical protocol for Straumann BLX implants, which are placed with the VeloDrill system (coronal widening). coDiagnostiX now uses the most efficient surgical protocol for all BLX implants.

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