

Release Notes 9.7.1

English

Feature updates

- 3DDX guided surgery kit: Support for new laterally opened sleeves added.
- Implant Solutions sleeve system updated.
- Neoss sleeve system updated to support 5mm and 6mm sleeves.
- Biologitech sleeve system support was added.
- Straumann guided surgery does now support rotationally locked sleeve holes.
- French, Japanese and Spanish translations updated.
- Updated the English user help.

Bug fixes

- Fixes problems with the digital drill guide wizard, dongles after computer standby, the screenshot printout, display of DICOM slices, and an online update problem with network dongle installations.

Deutsch

Neuerungen

- 3DDX Guided Surgery Kit: Unterstützung für neue, lateral geöffnete Hülse hinzugefügt.
- Implant Solutions' Hülsensystem wurde aktualisiert.
- Neoss' Hülsensystem wurde aktualisiert und unterstützt jetzt auch 5mm und 6mm Hülsen.
- Neues Hülsensystem: Biologitech
- Straumann-Guided-Surgery-Schablonen unterstützen nun rotationsgeschützte Hülsenlöcher.
- Die Übersetzung ins Französische, Japanische und Spanische wurde aktualisiert.
- Die Englische Benutzerhilfe wurde aktualisiert.

Fehlerbehebungen

- Behebt unter anderem Probleme im digitalen Bohrschablonenassistenten, im Screenshot-Ausdruck, mit Dongles und System-Standby, mit der Darstellung von DICOM-Daten und ein Online-Update-Problem im Zusammenhang mit Netzwerkdongle-Installationen.