

Release Notes 9.7.6

English

- Fixes a problem with plans that have been received via caseXchange. If such plans contained implants which were not available in the local implant database, some plans could not be loaded.
- Improved coDiagnostiX' ability to automatically fix problems in its system environment.

Deutsch

- Behebt ein Problem mit Plänen, welche über caseXchange empfangen wurden. Wenn diese Pläne Implantate enthielten, welche nicht in der lokalen Implantatdatenbank enthalten waren, konnten einige Pläne nicht geladen werden.
- Die coDiagnostiX-Funktion zum automatischen Beheben von Problemen an der Systemumgebung wurde verbessert.

Italiano

- Corretto il problema con le pianificazioni ricevute via caseXchange. Se le pianificazioni contenevano impianti che non erano disponibili nel data base locale degli impianti, alcune pianificazioni non si sarebbero potute caricare.
- Migliorate le capacità di coDiagnostiX di risolvere automaticamente i problemi nel suo ambito funzionale.